

Drytex® Playmaker®

Ideal Application

The Drytex Playmaker, suitable for mild to moderate ACL, MCL and LCL instabilities, is a neoprene alternative constructed of a unique nylon core and polyester Lycra® fabric that allows for improved airflow. The lightweight polyester-nylon weave fabric is ideal for patients who have neoprene-based irritations or allergies, or live in warm weather climates.

Features and Benefits

- Polycentric hinges with contoured condyle shells
- 4-point dynamic leverage system
- Extension stops: 10°, 20°, 30° and 40° (10° installed)
- Flexion stops: 45°, 60°, 75° and 90°
- Alternative to neoprene
- Lightweight, breathable fabric

Ordering Information

PART NUMBER

11-0557-x-06060

SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella

x=1=XS=13"-15 1/2" (33-39cm) x=2=S=15 1/2"-18 1/2" (39-47cm)

x=3=M=18 1/2"-21" (47-53cm) x=4=L=21"-23 1/2" (53-60cm)

x=5=XL=23 1/2"-26 1/2" (60-67cm) x=6=XXL=26 1/2"-29 1/2" (67-75cm)

Colors Available: black



NEVER STOP GETTING BETTER™

www.physiosupplies.com